Easy Power

Easy Power

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HEM – Hyperthermic Electro Massage

Technology based on the use of 2 hexapolar handheld applicators that act simultaneously by delivering different active energies that are synergistic.

Technically it is as if the operator were working with **two handpieces at the same time**.

The intensity adjustment of the various technologies can be managed through a remote control that the operator places on the wrist.

HEM is a multifunctional technology.

A highly technological treatment but at the same time perceived by the subject treated as a **warm cuddle, a massage where operator and technology interact in the application.**

It exploits the ability of very high frequency radio waves to produce heat in the tissues they pass through. The heat determines a contraction of the connective fibers which shorten and swell; as a result of this change, the fibroblasts are driven to produce new collagen and elastic fibers, as well as new hyaluronic acid and connective tissue substance. Neo fibrillogenesis is reflected in a significant increase in hydration, elasticity and tissue tone and is used to improve skin that has lost shape and firmness.

Secondly, the deep heat leads to a vasodilatation of the dermal layers, capable on the one hand of

reducing edema and swelling and on the other of oxygenating the tissues to the fullest. This ability is used to obtain a more nourished and luminous skin and to reduce water retention. In addition, the thermogenesis carried in depth, at the hypodermic level, determines an increase in local metabolism with an indirect lipolytic effect.

The endogenous heat generated also has a strong relaxing and decontracting effect on the muscles, thus optimally preparing the ground for the subsequent phases.


HEM also emits specific and studied electropulses: sudden pulsed electric fields capable of creating openings and channels through the epidermal cells. In this way any active ingredients applied are able to penetrate more easily in depth and perform their specific action better. This method, in the medical field, is also called **transdermal delivery**. In aesthetics, the applications are numerous and are linked to the active ingredient that is made to penetrate the skin; in general, antioxidant, lightening, draining, reducing, lipolytic etc. products are used.

DEMS – Diffuse Electro Myo Stimulation

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DEMS – Diffuse Electro Myo Stimulation

Method that uses electrical impulses of various frequency, intensity and duration to obtain a targeted contraction on certain muscle bundles, but also widespread on all minor, voluntary and involuntary muscle bundles present in the treated area. To obtain a correct and pleasant muscle movement for the treated subject, the waveform is biphasic compensated.

DEMS is applied to the muscles with an innovative system of electrically conductive bands which allow an ideal diffusion of stimulation in the treated area and allow the treatment to be concentrated on any areas with greater need. It spread the stimulation over multiple muscle bundles and make the physical activity generated more physiological. By varying the frequency, last of the pulse and the pause, the treatment is aimed more at aerobic and therefore slimming rather than anaerobic and therefore toning rather than microcirculatory activities. This technique allows you to perform a contraction without tiring (therefore for all ages) and make work both voluntary and involuntary fibers.

A DEMS phase carried out after a HEM Radio wave phase benefits from the endogenous heating generated in the area and allows to use high intensity electrical stimulation without discomfort.

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DES – Draining Electro Stimulation

*Technology essentially aimed at **venous and lymphatic drainage, extremely suitable for blemishes due to excessive stagnation of localized liquids**. It takes advantage of the electrical impulses to perform a peristaltic massage through the sequential activation of the stimulating channels, which create a significant increase in lymphatic migration towards the lymphatic drains.*

Essential treatment also to drain the body after phases aimed at slimming, toning and multifactorial treatment of cellulite.

*It is applied with an innovative system of **soaked electroconductive bands**.*

The treatment consists of two phases:

Phase A: characterized by electromagnetic stimulation capable of increasing peripheral circulation capillarization phase.

Phase B: characterized by a peristaltic massage performed through the sequential activation of the stimulating channels.

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Parti applicate: tipo BF / Applied parts: BF type

Guida operatore: display touch screen / Operator guide: touch screen display

Dimensioni: 57x48x50 cm / Overall dimensions: 57x48x50 cm

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